



Size and colour chart

Tights – Adult – Male Above Knee or Ankle

SIZES

Size	Waist (a)	Buttock (b)	Top leg (c)	Knee (d)	Inside leg length to ankle (e)	Inside leg length to knee (f)
XS	72 – 77cm	90 – 95cm	53 – 56cm	34 – 36cm	67.5cm	26cm
S	77 – 82cm	95 – 100cm	56 – 59cm	36 – 38cm	69cm	26.5cm
M	82 – 87cm	100 – 105cm	59 – 63cm	38 – 41cm	70.5cm	27cm
L	87 – 92cm	105 – 110cm	63 – 67cm	41 – 44cm	72cm	27.5cm
XL	92 – 97cm	110 – 115cm	67 – 71cm	44 – 47cm	73.5cm	28cm
XXL	97 – 102cm	115 – 120cm	71 – 75cm	47 – 50cm	75cm	28cm


FABRIC OPTIONS






MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest millimetre.
- Tape should be snug NOT tight.
- Waist (a) measurement is taken over the navel.
- Buttock (b) measurement is taken over the widest part of buttocks
- Knee (d) measurement is over the middle of the patella.
- Inside leg length to ankle (e) is measured to the end of the tights.
- Inside leg length to knee (f) is measured to above patella.



COLOURS



Please note garments *cannot* be altered or repaired by Second Skin.



Second Skin Pty Ltd
40 O'Malley Street
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455



Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937

secondskin.uk



RF ADULT-SC UK 2026-05



Size and colour chart

Tights – Adult – Female Above Knee or Ankle

SIZES

Size	Waist (a)	Buttock (b)	Top leg (c)	Knee (d)	Inside leg length to ankle (e)	Inside leg length to knee (f)
XS	67 – 72cm	92 – 97cm	54 – 57cm	34 – 36cm	67.5cm	26cm
S	72 – 77cm	97 – 102cm	57 – 60cm	36 – 38cm	69cm	26.5cm
M	77 – 83cm	102 – 108cm	60 – 64cm	38 – 41cm	70.5cm	27cm
L	83 – 89cm	108 – 114cm	64 – 68cm	41 – 44cm	72cm	27.5cm
XL	89 – 95cm	114 – 120cm	68 – 72cm	44 – 47cm	73.5cm	28cm
XXL	95 – 101cm	120 – 126cm	72 – 76cm	47 – 50cm	75cm	28cm


FABRIC OPTIONS






MEASUREMENTS

Measuring Tips:

- All circumferences to be measured to the nearest millimetre.
- Tape should be snug NOT tight.
- Waist (a) measurement is taken over the navel.
- Buttock (b) measurement is taken over the widest part of buttocks
- Knee (d) measurement is over the middle of the patella.
- Inside leg length to ankle (e) is measured to the end of the tights.
- Inside leg length to knee (f) is measured to above patella.



COLOURS



Please note that garments *cannot* be altered or repaired by Second Skin.



Second Skin Pty Ltd
40 O'Malley Street
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455



Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937

secondskin.uk





Second Skin

CUSTOMISED COMPRESSION GARMENTS & SPLINTS

Fabric options

All Second Skin garments are made from a range of medical grade stretch fabrics. Several highly durable medical grade fabrics are incorporated in our garment designs, to address the individual needs of each client.

POWERNET FABRICS

Powernet is an open weave fabric which is prescribed for healed skin and for areas where additional support and compression is required. This fabric is particularly durable in work, school or outdoor activities. Powernet is often used in burns scar management and lymphoedema.

Our Powernet fabric now incorporates patented cooling properties, HeiQ Smart Temp, which provides outstanding thermal comfort to clients by keeping them cooler and dryer while wearing their garment.

Powernet fabric is available in light tan, dark tan and black.

SHIMMER FABRICS

Shimmer is a closed weave fabric which provides a softer interface with the skin. Shimmer fabric contours well and is comfortable and effective over joint surfaces and areas of skin fragility such as the elbow, ankle and knee. Shimmer fabric is used for early burns and flap management, external fixators, vascular conditions and certain types of lymphoedema.

Shimmer fabric is available in light tan, dark tan and black.

HYDROPHOBIC FABRICS

Hydrophobic fabric is used as both a lining and a base fabric in compression garments. This fabric wicks moisture away from the surface of the skin to evaporate, protecting the skin surface and ensuring optimal comfort.

Our hydrophobic fabric has an antimicrobial component embedded within the matrix of the fabric. This ensures our products are breathable and highly effective at removing bacteria from the fabric surface, to provide the wearer with all year round comfort. Due to this innovation, our garments are frequently prescribed to assist with the management of repeated skin breakdown and fragility.

Hydrophobic fabric is available in light tan, dark tan and black.

PERTH OFFICE

40 O'Malley Street
OSBORNE PARK
WA 6017

P: +61 8 9201 9455

E: perth@secondskin.com.au

UK OFFICE

Unit 10, Spectrum House
32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP

P: 0203 6279937

E: london@secondskin.com.au

 Follow us on LinkedIn