

Size and colour chart

Knee-High Sock – Child Open or Closed Toe


SIZES

Size	Prox end circumference (a)	Dorsal ankle circumference (b)	Foot length (c)	Height from floor (d)
2	21.1 – 22.9cm	15.6 – 18.4cm	13.1 – 14.9cm	18.5cm
4	22.9 – 24.7cm	18.4 – 21.2cm	14.9 – 16.7cm	21.5cm
6	24.7 – 26.4cm	21.2 – 24cm	16.7 – 18.5cm	25.0cm
8	26.4 – 28.1cm	24 – 26.8cm	18.5 – 20.3cm	28.5cm
10	28.1 – 29.9cm	26.8 – 30.1cm	20.3 – 22.6cm	33.0cm
12	29.9 – 31.6cm	30.1 – 33.2cm	22.6 – 24.7cm	36.5cm


MEASUREMENTS

Measuring Tips:




- All circumferences and lengths to be measured to the nearest millimetre.
- Tape should be snug NOT tight



COLOURS



FABRIC OPTIONS



Please note that garments *cannot* be altered or repaired by Second Skin.



Second Skin Pty Ltd
40 O'Malley Street
Osborne Park WA Australia 6017
Email: perth@secondskin.com.au
Phone: +61 8 9201 9455



Second Skin
10/32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP
Email: london@secondskin.com.au
Phone: +44 203 627 9937

secondskin.uk





Second Skin

CUSTOMISED COMPRESSION GARMENTS & SPLINTS

Fabric options

All Second Skin garments are made from a range of medical grade stretch fabrics. Several highly durable medical grade fabrics are incorporated in our garment designs, to address the individual needs of each client.

POWERNET FABRICS

Powernet is an open weave fabric which is prescribed for healed skin and for areas where additional support and compression is required. This fabric is particularly durable in work, school or outdoor activities. Powernet is often used in burns scar management and lymphoedema.

Our Powernet fabric now incorporates patented cooling properties, HeiQ Smart Temp, which provides outstanding thermal comfort to clients by keeping them cooler and dryer while wearing their garment.

Powernet fabric is available in light tan, dark tan and black.

SHIMMER FABRICS

Shimmer is a closed weave fabric which provides a softer interface with the skin. Shimmer fabric contours well and is comfortable and effective over joint surfaces and areas of skin fragility such as the elbow, ankle and knee. Shimmer fabric is used for early burns and flap management, external fixators, vascular conditions and certain types of lymphoedema.

Shimmer fabric is available in light tan, dark tan and black.

HYDROPHOBIC FABRICS

Hydrophobic fabric is used as both a lining and a base fabric in compression garments. This fabric wicks moisture away from the surface of the skin to evaporate, protecting the skin surface and ensuring optimal comfort.

Our hydrophobic fabric has an antimicrobial component embedded within the matrix of the fabric. This ensures our products are breathable and highly effective at removing bacteria from the fabric surface, to provide the wearer with all year round comfort. Due to this innovation, our garments are frequently prescribed to assist with the management of repeated skin breakdown and fragility.

Hydrophobic fabric is available in light tan, dark tan and black.

PERTH OFFICE

40 O'Malley Street
OSBORNE PARK
WA 6017

P: +61 8 9201 9455

E: perth@secondskin.com.au

UK OFFICE

Unit 10, Spectrum House
32-34 Gordon House Road
Gospel Oak, LONDON NW5 1LP

P: 0203 6279937

E: london@secondskin.com.au

 Follow us on LinkedIn